



18	1	Skydd för öppning	551112
17	1	Dockningsstation	551101
16	2	Löpbana GH1	551099
15	1	Laddkabel GH1	551100
14	1	Transformator	550200
13	1	Skyddskåpa, botten	551102
12	1	Ni-mH Batteri pack 24V/2.8Ah	550574
11	1	Lyftband	550547
10	2	Koppling löpbana GH1	551103
9	1	Manöverdosa GH1	551104
8	1	Bandgenomföring inkl switch	550564
7	1	Motor	551105
6	1	Drivrem	550580
5	2	Etikett botten GH1	551106
4	1	Laddkort PCB	550573
3	1	Huvud kretskort	550548
2	1	Skyddskåpa, topp	551107
1	2	Sidoskydd, set 2 + 2 st	551108
Pos no.	Antal	Artikel	Ritning/Artnr.

Tolerance:		Surface treatment:		
Material:		Colour:	Weight (kg):	
Designed by: ob	Drawing by: ob	Date: 01-12-2010	Scale: 1:5	Sheet: A3
<b>Guldmann</b> <small>Confidential: Property of V. Guldmann A/S. Not to be handed over to, copied or used by third party.</small>		V. Guldmann A/S Graham Bells Vej 21-23A DK-8200 Århus N Tel. +45 8741 3151 Fax +45 8741 3131 E-mail: info@guldmann.com WWW.guldmann.com		Change no.:
		Product: GH1 Part: GH1 Spare parts Drawing no.: 55225X		Edition: 0