



16.1	1	553654	Batteribox, GL5, 155kg
16.2	1	553653	Batteribox, GL5, 205kg
16.3	1	553854	Kretskort, GL5, 205kg
16.4	1	553853	Kretskort, GL5, 155kg
16.5	1	554012	Batteri, kit bestående av 4 st GL5
16.6	1	554011	Switch, Nödstopp GL5
16.7	1	554010	Switch, Nödsänkning GL5
15	1	553860	Bult för sockel, mast GL5
14	1	553859	Underrede, GL5
13	1	553858	Manöverdosa, GL5
12	1	553857	Skydd, gummi 30
11	1	553856	Elkabel laddare, Grå GL5
10	1	553855	Laddare GL5
9	1	553852	Körhandtag, GL5
8	1	553851	Breddningsmekanism GL5
7.2	1	553850	Hjul, 80mm
7.1	1	553849	Hjul, 100mm
6.2	1	553848	Hjul, broms, 80mm
6.1	1	553847	Hjul, broms, 100mm
5	1	553846	Infästningsdetaljer ben/sockel GL5
4	1	553845	Ben GL5
3	1	553844	Actuator, benbreddning
2.2	1	553843	Actuator, lyft
2.1	1	553842	Actuator, lyft, manuell nödsänkning
1.2	1	553841	Mast & Lyft arm
1.1	1	553840	Mast & Lyft arm, Manuell nödsänkning
Pos nr.	Antal	Artnr	Benämning

Tolerance:		Surface treatment:		
Material:		Colour:	Weight Kg :	
Designed by: jta	Drawing by: jta	Date: 17-05-2011	Scale:	Sheet: A3
<b>Guldmann</b> <small>Confidential: Property of V. Guldmann A.S. Not to be handed over to, copied or used by third party.</small>		V. Guldmann A S Graham Bells Vej 21-23A DK-8200 crrhus N Tel. +45 8741 3151 Fax +45 8741 3131 E-mail: info@guldmann.com WWW.guldmann.com		Change no.:
		Product: Mobile Hoist Part: Spareparts - GL5		Drawing no.: Edition: 1